


MODULAR, FLEXIBLE, UPGRADEABLE

ROBOTIC SOLUTIONS





THE ULTIMATE DESIGN

The Synchron Xperimate is a truly modular state-of-the-art robotic system, combining sophisticated and flexible liquid handling with robotic manipulation. The modular design allows the hardware to be configured to suit a very wide range of applications.


The modular Xperimate is a perfect platform for all applications, and simplifies both system integration and stand alone robotic automation. The sleek, uncomplicated design integrates robotic and liquid handling functionality into the robotic arm, which glides on an X-rail to access the large flexible deck, and integrated or surrounding modules.

STAND-ALONE WORKSTATION

Your needs for automation may change during the long lifetime of your Xperimate robot. The unique modular design means that you can select suitable elements and performance characteristics to make up the ideal workstation, and then upgrade later to increase throughput or to more fully automate a protocol.

FULL INTEGRATION


With all functionality housed in the arm, this entire module can be easily integrated into an instrument project. Our CAN bus protocol simplifies integration. Alternatively, the Xperimate workstation, including deck, can be customized and optimized as a private label platform to automate any pipetting protocol.



Xperimate is composed of a basic frame with one or more robotic arms. Each arm is constructed from one or two arm sections, and each section can carry 1, 2 or 4 pipetting tips, or a robotic handling tool. Several standard configurations are available, and Xperimate can be customized for all applications.



CUSTOMIZING CONFIGURATIONS



By selecting from a wide range of robotic modules, functional "Plug and Play" modules, options and accessories, you can customize your Xperimate to suit the evolving needs of your laboratory. And in most cases you can change the configuration in the future by upgrading.

DESIGN YOUR OWN PIPETTING ARM

- Up to 8 tips per arm
- Flexispan = individual Y and Z movement
- · Individual level sensors
- Interchangeable washable and disposable tips.

ROBOTIC HANDLING


- Separate arm, same X-slide
- · Software controlled for safety and efficiency
- 270° rotating microplate gripper
 - turn a plate for cross pipetting
 - improved access to other modules on or off the deck
- Increase deck capacity by stacking plates
- Increase efficiency by moving plates between functional modules
 - incubators, coolers, shakers, readers and washers, thermocycler, vacuum assembly, centrifuge etc.

DECK LAYOUT


- 3 standard X-ranges, and 3 standard deck widths (approx. 100 cm, 150 cm and 200 cm)
- Extended X-range option
 - reach modules outside deck
- $\bullet \ \ \text{Decktrays hold samples, plates, stackers, shakers, reagents and the tip wash station}$
- Simple to use drag and drop software helps you define the labware layout.







Stacking plates



Decktrays

STANDARD XPERIMATE CONFIGURATIONS

	ARM 1		ARM 2		
Model	Section 1	Section 2	Section 1	Section 2	
40 00	4 pipetting tips	-	_	_	
40 0H	4 pipetting tips	_	-	robotic handler	
44 00	8 pipetting tips		-	_	
44 0H	8 pipetting tips		-	robotic handler	
44 44	8 pipet	8 pipetting tips		8 pipetting tips	

The above table shows some examples of Xperimate configurations. Xperimate units are available in various deck widths and with optional X-range extension. Other models may be available.



PRINCIPLES OF LIQUID HANDLING

The robotic pipetting arm module encloses everything necessary to perform the liquid handling, reducing the need for space consuming external pumps and extended lengths of feed tubing. Xperimate uses a principle of positive displacement of liquid in a fluid filled system, where the aliquots of sample and reagent are automatically separated by an airgap. Xperimate can pipette using all or only some of the individual channels together, either in single aliquots or by multi-dispensing.



- · Mini toothed wheel pump
- Pipettes 0.5 ul to 2.3 ml, plus continuous system liquid dispensing
- · Excellent precision
- No maintenance, no replacement parts.

PIPETTING TIPS

- Up to 8 pipetting tips per arm, up to 2 arms per standard Xperimate
- Interchangeable washable and disposable tips
 - select tip format and size for each protocol step
- Flexispan independent YZ drive and level detector for each tip
 - access consecutive or unevenly spaced tubes and wells (e.g. for cherry picking)
- Each tip can exert a force of 30 N (3 kg)
 - ideal for piercing septa.



- · Unique design optimized for precision and minimal carry-over
 - various coatings available
- · Total pipetting volume 2.3 ml
 - sample or reagent contained in the tip/tubing
 - smooth interfaces optimizes precision.




- · For applications requiring zero carryover
- Various formats and volumes.



Micro-pump

Xperimate arm assembly



Flexispan



Individual tip level detectors



Washable pipetting tip




Disposable tips




ADDING A THIRD DIMENSION

The Xperimate deck is large and open on all sides. It provides a canvas on which to place labware and optional modules.

Sample tube racks, microplates and reagents are loaded on the slimline decktrays, which accurately position the labware ready for automation. Most Xperimate accessories and options are simply placed on the deck and can be removed when not in use, or rearranged to suit another application protocol.







ROBOTIC HANDLING

The separation of the robotic and pipetting arms ensures that robotic function can work independently to increase throughput and efficiency. Plates and other labware are quickly and safely carried between stackers, incubators and other functional modules. Plates can be turned for efficient cross pipetting tasks. When installed on a Xperimate with extended X-range, the robotic arm can access modules and accessories situated outside the Xperimate deck.

- · Creates a useful third dimension to the deck
- Transport plates with or without their lids; remove and replace lids
- Pipette to a plate on the top of a stack
- 270° rotating microplate gripper.





The narrow microplate gripper can access the loading bay of the centrifuge.



EXPANDING THE RANGE OF AUTOMATION

The annual and a second and a s

The basic ID-trax module holds 96 samples and can be expanded in batches of 96.

The functionality of the Xperimate can be expanded by using specialized modules and accessories. Most Synchron modules are simply placed on the deck and are fully Plug and Play. Third party modules may also be integrated.

SAMPLE TRACKING - ID-TRAX


- · Increase security with real-time ID tracking
- · Reads barcodes on sample tubes, plates and reagents
- Selective (worklist) pipetting, or batch mode
- Decktop Plug and Play module can be removed to increase deck capacity when not in use.



Incubator

INCUBATOR

- Holds up to 4 plates, with or without covers
- · Plug and Play
- · Individual plate drawers
- · Minimize drafts and edge effect
- Incubation time and temperature are software controlled.



ixion

IXION - INTEGRATED MICROPLATE CENTRIFUGE

Accommodating 2 microplates, deep well plates or gel cards, the lxion can be integrated into the deck of the Xperimate, or can be positioned alongside. Loading, unloading and functions of the lxion are controlled by the Xperimate software.

Features


- rcf up to 2,000 g
- Centrifugation time 5 seconds to 9 hours
- · Loading shutter opens and closes under software control
- Automated homing position to permit robotic loading and unloading
- Motion is prevented if the lid is not fully shut
- Automatic sensor detects imbalance, and stops centrifuge if there is excess vibration.

OTHER MODULES

The Xperimate design utilizes a CAN protocol to facilitate integration of third party products. The deck is large, and has few limitations for physical integration of hardware. Additionally the Xperimate can be delivered with an extended X-range, so that the robotic arms can reach outside the deck to access modules placed alongside.




X-AP – SOFTWARE FOR EVERY APPLICATION



Deckmate allows you to easily define the Xperimate deck layout



X-AP stores your pre-optimized methods and displays them in a list for easy selection



X-AP is powerful and very flexible software designed to set-up and automate all your liquid handling tasks. It can be used to communicate with and/or control functional modules for complete assay automation.

DECKMATE

A drag and drop software used to design the arrangement of samples, reagents, plates, etc. on the deck for each protocol. Deckmate is delivered with an extendable library of labware, so you can focus on designing the best deck for your application.

X-AP

A graphic user interface used to compose and optimize protocols. X-AP has multiple user levels. The programming level allows the definition and optimization of every step. Simply transpose your application step by step into a X-AP method. To automate complex tasks, you can combine several methods into a Script.

- * X-AP utilizes a robust database concept which keeps track of all data, and provides a seamless interface with Microsoft $^{\circ}$ Excel
- X-AP includes a service utilities program, X-Util, designed to simplify system diagnostics.
 X-Util can also be used to fully customize unique applications, by addressing each module and function individually
- A powerful simulator software is also available compose and test your methods prior to running them on the Xperimate.

Scheduling


- · Online scheduling constantly re-evaluates timelines during a run
- · Actions shared by different protocols can be consolidated
- · Add or remove actions while an application is running
- Time calculation automatic timing of actions, and correction of estimates
- · Graphic display
- Active-X control can be used for external modules.

Connectivity

- DDE data exchange
- Integrate external module Dlls and "call" from within a run
- Additional functions can be executed via the optional I/O board.
- Import/Export from Excel methods, labware and pipetting lists can be imported or exported before, during or after a run.

APPLICATION SOFTWARE

X-AP is a very flexible software that can be used to precisely define and optimize your applications. Some pre-optimized application packages are also available (e.g. ELISA, blood grouping, IFA etc.). Please contact us for more details.



The optional ELISA user interface simplifies the complete automation of plate preparation, processing and data reduction.



SPECIFICATIONS

ACCESSO	DRIES
---------	-------

Microplate shaker Capacity 1 MTP or deep-well plate; linear, orbital and customized shaking modes

Decktrays A decktray holds 96 sample tubes; or up to 5 microplates/stacks of 10; DT racks; or reagent racks etc

The deck model 200 can accommodate 13 trays The deck model 150 can accommodate 9 trays The deck model 100 can accommodate 5 trays

(Tip access to some wells of the 4 & 5 plate trays may be limited)

PIPETTING PERFORMANCE

Carry over Washable tip: down to 10⁻⁷ with sufficient tip wash; Disposable tip: Zero

Liquid/clot detector Capacitance

Precision – gravimetric evaluation Vol (ul) CV (optimized conditions)

10 1 % > 20 < 0.5 %

Pumps Individual mini-toothed wheel pump for each tip

Tip wash

Tip wash station can be placed anywhere on the deck

Speed/throughput Depends on Xperimate configuration - please consult Synchron

Tip design Washable: Optimized for precision and minimal carry-over, various coatings available

Disposable: 10 ul, 100 ul, 1000 ul, various formats

Minimum equal spacing between tips 9 mm Maximum equal spacing of tips 55 mm

Maximum single spacing 350 mm

Volumes $0.5~\mu l - 2300~\mu l$. Volumes outside this range may be possible

ROBOTIC PERFORMANCE

Flexispan

Force in z-direction 30 I

Arm force 240 N at the far end Handler Can lift 3 kg

Position Accuracy +/- 0.2 mm (384 well plates)

z-travel range 150 mm

ID-TRAX (barcode reader)

CCD camera reads most commercially available barcodes, max digits 20

Standard capacity 96 samples, extendable in batches of 96

INCUBATOR TOWER

Capacity 4 microplates. Temperature controllable from 5° C above ambient temperature

to 49° C, in steps of 0.1° C, accuracy 0.2° C

OPERATING CONDITIONS

Temperature 15° C to 32° C Relative humidity 10-85% at 32° C

ELECTRICAL AND COMMUNICATION Plug and Play modules simplify the composition of an application system

RS232, multiple I/O ports, and a CAN bus facilitate systems integration

Power Self regulating power supply, 85 to 264 volts AC; 350 watts

SAFETY AND SECURITY

Safety Operator safety is ensured by the safety guard or an optional cover

Xperimate is designed for IEC 1010 compliance and bears the CE mark

Security Sophisticated software combined with real time barcode tracking can eliminate sample tracking errors

Software user level definition prevents unauthorized changing of performance data

Processes and actions are monitored and logged



ABOUT SYNCHRON

Synchron was founded in 1985, based on many years of experience developing and marketing flexible robotic liquid handling and automation. Synchron offers expertise in many areas of automation for clinical, blood bank, pharmaceutical and research applications. Our broad understanding of the need for more flexible systems and modules is applied to the development of products for the most demanding laboratory applications. Synchron is a dedicated and experienced partner who can help you realize your automation projects efficiently and cost effectively.

